

# ALK5 inhibitor I BI-4659

BI-4659



# **Table of contents**

Summary	2
Chemical Structure	
Highlights	3
Target information	3
In vitro activity	4
In vitro DMPK and CMC parameters	4
In vivo DMPK parameters	
Negative control	5
Selectivity	6
Co-crystal structure of the Boehringer Ingelheim probe compound and the target protein	6
Reference molecule(s)	6
Supplementary data	7
References	8



#### **Summary**

BI-4659 inhibits ALK5 with an IC $_{50}$  value of 19 nM and blocks the cellular phosphorylation of SMAD2/SMAD3 with an EC $_{50}$  of 185 nM. BI-4659 is a suitable tool for testing biological hypotheses *in vitro*.

#### **Chemical Structure**

Figure 1: 2D structure of BI-4659, an ALK5 inhibitor




Figure 2: BI-4659, 3D conformation

#### **Highlights**

BI-4659 is a potent ALK5 (TGF $\beta$ R1) inhibitor (IC $_{50}$  = 19 nM). This compound showed good selectivity against a broad panel of other kinases. BI-4659 is able to block the phosphorylation of SMAD2 and SMAD3 in cells and is a suitable tool for in vitro studies $^1$ . However, it is not suited for in vivo experiments due to the very high clearance abserved in rats.

#### **Target information**

Transforming growth factor  $\beta$  (TGF $\beta$ ) is a pluripotent cytokine involved in the regulation of various biological processes such as cell proliferation, differentiation, migration, adhesion, apoptosis, and epithelial-to-mesenchymal transition (EMT). Therapeutic approaches to inhibit its signaling by targeting TGF $\beta$  receptor I (TGF $\beta$ RI/ALK5) are discussed for the treatment of diseases such as Idiopathic Pulmonary Fibrosis (IPF) and cancer<sup>2,3,4</sup>.




Figure 3: Crystal structure of ALK5 complexed with a close analog of BI-4659 (PDB code 2X7O)<sup>1</sup>

# In vitro activity

PROBE NAME / NEGATIVE CONTROL	BI-4659	BI-4101
MW [Da, free base] <sup>a</sup>	440.5	418.5
Inhibition of ALK5 (Kinase Glow assay, IC <sub>50</sub> ) [nM] <sup>b</sup>	19	>50000
Inhibition of phosphorylation of SMAD2 and SMA3 in HaCaT cells (EC50) [nM]	185	n.d.

<sup>&</sup>lt;sup>a</sup> For the salt form you will get, please refer to the label on the vial and for the molecular weight of the salt, please refer to the FAQs

# *In vitro* DMPK and CMC parameters

PROBE NAME / NEGATIVE CONTROL	BI-465	59		BI-4101
logD @ pH 2/ 11	0.7 / 3	.1		n.d.
Solubility @ pH 6.8 [μg/mL]	2			n.d.
Solubility @ pH 4 [µg/mL]	84			n.d.
Solubility @ pH 2 [µg/mL]	71			n.d.
Caco-2 permeability AB @ pH 7.4 [*10 <sup>-6</sup> cm/s]	8			n.d.
Caco-2 efflux ratio	6			n.d.
Microsomal stability (human/mouse/rat) [% QH]	26	23	<22	n.d.
hERG [inh. % @ 10 μM]	38.5	•		n.d.
hERG IC <sub>50</sub> [μM]	>10			n.d.

 $<sup>^{\</sup>mathrm{b}}$  For detailed assay conditions see reference 1

CYP 3A4 (IC <sub>50</sub> ) [μM]	>50	n.d.
CYP 2C8 (IC <sub>50</sub> ) [μM]	>50	n.d.
CYP 2C9 (IC <sub>50</sub> ) [μM]	>50	n.d.
CYP 2C19 (IC <sub>50</sub> ) [μM]	>50	n.d.
CYP 2D6 (IC <sub>50</sub> ) [μM]	>50	n.d.

#### In vivo DMPK parameters

BI-4659 showed very high clearance in rat and therefore is not suited if in vivo studies

BI-4659	RAT
Clearance [%Q <sub>H</sub> ] <sup>a</sup>	123
Mean residence time after i.v. dose [h] <sup>a</sup>	1.7
t <sub>max</sub> [h] <sup>b</sup>	3.3
C <sub>max</sub> [nM] <sup>b</sup>	27.7
F [%] <sup>b</sup>	50
Vss [L/kg] <sup>a</sup>	9

<sup>&</sup>lt;sup>a</sup> *i.v. dose:* 4.5 mg/kg <sup>b</sup> *p.o.* dose: 45 mg/kg

### **Negative control**

BI-4101 shows no inhibition of ALK5 in the Kinase Glow assay (IC<sub>50</sub> > 50  $\mu$ M) and therefore is a suitable negative control for *in vitro* experiments (BI-4101 = Compound **48**<sup>1</sup>)

Figure 4: BI-4101 which serves as a negative control

#### Selectivity

BI-4659 shows good selectivity in kinase panels. The compound shows no inhibition of 218/232 kinases tested at 2  $\mu$ M. Cross-reactive kinases (%inhibition data at 0.2  $\mu$ M): ABL1 (67), BLK (85), CSF1R (FMS) (82), FGR (91), FLT3 (66), FYN (74), LCK (100), LYN A (79), LYN (78), MAP4K5 (KHS1) (77), MELK (80), NTRK3 (TRKC) (78), RET (79), SNF1LK2 (96), YES1 (90).

Selectivity data is available on opnMe.com for download free of charge.

SELECTIVITY DATA AVAILABLE	BI-4659	BI-4101
SafetyScreen44™ with kind support of	Yes	Yes
Invitrogen®	Yes	No
DiscoverX®	Yes	No
Dundee	Yes	No

# Co-crystal structure of the Boehringer Ingelheim probe compound and the target protein

No in-house structure is available for BI-4659, but for related compound (PDB code 2x7o).

# Reference molecule(s)

SB-505124, SB-525334, GK6604, SD-208, LY-2157299, EW-7197, GW788388 and others. For a review on ALK5 kinase inhibitors in oncology $^5$ .

# Supplementary data

Selectivity data and 2D structure files can be downloaded free of charge from opnMe.



#### References

- Roth G. J., Heckel A., Brandl T., Grauert M., Hoerer S., Kley J., Schnapp G., Baum P., Mennerich D., Schnapp A. and Park J. E. Design, Synthesis, and Evaluation of Indolinones as Inhibitors of the Transforming Growth Factor β Receptor I (TGFβRI) *J. Med. Chem.* 2010, 53, 7287-7295. DOI: 10.1021/jm100812a, PubMed.
- 2. De Gouvill A-C., Hue S. Inhibition of ALK5 as a new approach to treat liver fibrotic diseases *Drug News Perspect* **2006**, *19*, 85-90. DOI: 10.1358/dnp.2006.19.2.977444, PubMed.
- 3. Fuchs O. Inhibition of TGF- Signaling for the Treatment of Tumor Metastasis and Fibrotic Diseases *Current Signal Transduction Therapy* **2011**, *6*, 29-43. DOI:10.2174/157436211794109398.
- 4. Derynck R., Zhang Y. E. Smad-dependent and Smad-independent pathways in TGF-beta family signaling *Nature* **2003**, *425*, 577-584. DOI: 10.1038/nature02006, PubMed.
- 5. Ling L. E., Lee W-C. TGF-Beta Type I Receptor (ALK5) Kinase Inhibitors in Oncology *Curr. Pharm. Biotechnol.* **2011**, *12*, 2190-2202. DOI: 10.2174/138920111798808257, PubMed.
- Baum P., Schmid R., Ittrich C., Rust W., Fundel-Clemens K., Siewert S., Baur M., Mara L., Greunenbaum L., Heckel A., Eils R., Kontermann R. E., Roth G. J., Gantner F., Schnapp A. Park J. E., Weith A., Quast, K., Mennerich D. A Strategy to Qualify Chemical Compounds during Hit-to-Lead and/or Lead Optimization *PLoS One* 2010, 5, e14272. DOI: 10.1371/journal.pone.0014272, PubMed.

