

Complex I - ROS modulator

BI-4500

Table of contents

Summary	2
Chemical Structure.....	2
Highlights.....	3
Target information.....	3
<i>In vitro</i> activity.....	4
<i>In vitro</i> DMPK and CMC parameters	5
<i>In vivo</i> DMPK parameters.....	6
Negative control.....	6
Selectivity.....	7
Supplementary data	7
References.....	7

Summary

BI-4500 is a potent and selective NADH:ubiquinone oxidoreductase (Complex I) reactive oxygen species (ROS) modulator. BI-5330 is available as negative control.

Chemical Structure




Figure 1: 2D structure of BI-4500




Figure 2: Low-energy 3D conformation of BI-4500

Highlights

BI-4500 is a potent and selective NADH:ubiquinone oxidoreductase (Complex I) reactive oxygen species (ROS) modulator that reduces oxidative stress by inhibiting mitochondrial ROS generation. It preserves NADH:ubiquinone oxidoreductase activity, distinguishing it from classical Complex I inhibitors.

Target information

Selective inhibitors of mitochondrial Complex I–derived reactive oxygen species (ROS) represent an advanced strategy to control oxidative stress without disrupting essential cellular energy production¹⁻⁴. Mitochondrial Complex I (NADH:ubiquinone oxidoreductase) is a key component of the electron transport chain and plays a critical role in ATP synthesis. At the same time, it is a major intracellular source of ROS, particularly under conditions of stress or dysfunction, where electron leakage can lead to excessive superoxide formation.

Unlike traditional Complex I inhibitors, which suppress enzymatic activity and impair oxidative phosphorylation, selective ROS modulators act more precisely. They target specific sites within Complex I that are responsible for pathological ROS generation, reducing electron leakage while preserving normal electron flow and ATP production⁵⁻⁶. This selective mechanism allows for the attenuation of oxidative damage without compromising cellular metabolism.

Such compounds are powerful research tools to investigate molecular mechanisms underlying diseases characterized by mitochondrial dysfunction and elevated oxidative stress, including neurodegenerative disorders like Parkinson's and Alzheimer's disease, as well as cardiovascular and inflammatory conditions²⁻⁵. By lowering ROS levels at their source, these inhibitors may help protect cellular components such as lipids, proteins, and DNA from oxidative damage, thereby preserving tissue integrity and function. In addition, many of these molecules are designed with favorable pharmacokinetic properties, including sufficient stability, tissue distribution, and in some cases, the ability to cross the blood–brain barrier. This brain penetrance is particularly important for targeting central nervous system disorders. Overall, selective inhibitors of Complex I–derived ROS represent valuable research reagents, combining high specificity with a potentially improved safety profile compared to non-specific antioxidants or conventional mitochondrial inhibitors.



Figure 3: Cryo-EM structure of the multi-subunit Complex-I from human mitochondria (PDB code: 9I4I). The binding location of BI-4500 is unknown.

***In vitro* activity**

BI-4500 inhibits purified bovine Complex I mediated ROS-formation with an IC_{50} of 80 nM. Compounds from this structural class have been shown to have no species selectivity. The Complex I ROS modulator was also shown to be selective versus other ROS-generating enzymes (e.g. xanthine oxidase, monoamine oxidase and NO-synthases all $>10\mu M$) and does not show cytotoxicity at 30 μM in neuroblastoma cells, indicating that mitochondrial respiration is unaffected (modulation concept). Furthermore, the small molecule shows protection of glutathione-depleted neuroblastoma cells in a model of oxidative stress (HT22 cells, EC_{50} = 80 nM) and dopaminergic neurons in the MPTP-model, a neurotoxin inhibiting mitochondrial Complex I resulting in diminished ubiquinone-reduction and ATP synthesis as well as stimulated ROS-formation.

Probe name / Negative control	BI-4500	BI-5330
MW [Da] ^a	357.45	283.3
Complex I (IC ₅₀) (bovine) [nM] ^b	80	>10,000
HT22 cell protection (EC ₅₀) [nM] ^b	80	>10,000
Selectivity factor versus Xanthinoxidase (ROS generation)	>100	-

^a The molecule is supplied in salt form; for the molecular weight of the salt, please refer to the vial label.

^b see US2021061761

In vitro DMPK and CMC parameters

BI-4500 has acceptable solubility in water at neutral pH and high permeability in Caco2 and MDCK assays.

Probe name / Negative control	BI-4500	BI-5330
logD @ pH 7.4	2.4	0.3
Solubility @ pH 6.8 [µg/mL]	100	8
Caco-2 permeability AB @ pH 7.4 [$\times 10^{-6}$ cm/s]	69.3	<16
Caco-2 efflux ratio	0.9	-
MDCK permeability PappAB @ 1µM [10^{-6} cm/s]	65	-
MDCK efflux ratio	0.9	-
Microsomal stability (human/mouse/rat) [% Q _H]	<23 / 34 / <22	<66 / 76 / 38
Hepatocyte stability (human/mouse/rat) [% Q _H]	18 / - / 41	- / - / 72
Plasma Protein Binding (human/mouse/rat) [%]	80 / - / 81	-
hERG IC ₅₀ [µM]	> 10	-
CYP 3A4 (IC ₅₀) [µM]	12	>50
CYP 2C8 (IC ₅₀) [µM]	11	>50
CYP 2C9 (IC ₅₀) [µM]	>50	>50
CYP 2C19 (IC ₅₀) [µM]	46	-
CYP 2D6 (IC ₅₀) [µM]	19	>50

In vivo DMPK parameters

PK properties in several animal species are suitable for once or twice daily oral dosing of BI-4500 in acute or sub-chronic *in vivo* experiments. Initial safety studies in rats showed no cardio-vascular effects up to a 50 mg/kg dose.

BI-4500	Mouse ^a	Rat ^b
Clearance [% Q _H] ^a	14	27
Mean residence time after <i>i.v.</i> dose [h] ^a	3.3	3.3
t _{max} [h] ^b	1	0.8
C _{max} [nM] ^b	1,420	4,452
F [%] ^b	100	100
V _{ss} [l/kg] ^a	2.4	3.6

^a *i.v.* dose 3.6 mg/kg, *p.o.* dose 1.5mg/kg

^b *i.v.* dose 3.6mg/kg, *p.o.* dose 3.6 mg/kg

Negative control

BI-5330 is a structural analog and can be used as a negative control.




Figure 4: BI-5330 which serves as a negative control

Selectivity

BI-4500 showed high selectivity at 10µM versus a panel of receptors. It inhibited 4 out of 85 targets tested in SafetyScreen44™ by ≥ 50% (Free Radical Scavenging, 5-LO, L-peroxidase, 15-LO). 7 targets have an IC₅₀ <1µM in Invitrogen assay (SGKL, PLK3, MAPKAPK3, KDR, IKBKB, FGR, CAMK2A).

Neg. control BI-5330 showed high selectivity at 10µM versus a panel of receptors. It inhibited only 1 out of 43 targets in SafetyScreen44™ by >50% (PDE4D2).

Selectivity data available	BI-4500	BI-5330
SafetyScreen44™ with kind support of  eurofins	Yes	Yes
Invitrogen	Yes	No

Supplementary data

2D structure files can be downloaded free of charge from [openMe](https://openme.boehringer-ingenheim.com).

References

1. Murphy M. P. How mitochondria produce reactive oxygen species *Biochem J* **2009**, 417(1), 1–13. [DOI: 10.1042/BJ20081386](https://doi.org/10.1042/BJ20081386), [PubMed: 19061483](https://pubmed.ncbi.nlm.nih.gov/19061483/).
2. Janssen R. J. R. J., Nijtmans L. G., van den Heuvel L. P., Smeitink J. A. M. Mitochondrial complex I: Structure, function and pathology *J Inherit Metab Dis* **2006**, 29(4), 499–515. [DOI: 10.1007/s10545-006-0362-4](https://doi.org/10.1007/s10545-006-0362-4), [PubMed: 16838076](https://pubmed.ncbi.nlm.nih.gov/16838076/).
3. Koilkonda R. D., Guy J. Leber's Hereditary Optic Neuropathy-Gene Therapy: From Benchtop to Bedside *J Ophthalmol* **2011**, 2011, 179412. [DOI: 10.1155/2011/179412](https://doi.org/10.1155/2011/179412), [PubMed: 21253496](https://pubmed.ncbi.nlm.nih.gov/21253496/).
4. Piccoli C., Ripoli M., Quarato G., Scrima R., D'Aprile A., Boffoli D., Margaglione M., Criscuolo C., Michele G. de, Sardanelli A., Papa S., Capitanio N. Coexistence of mutations in PINK1 and mitochondrial DNA in early onset parkinsonism *J Med Genet* **2008**, 45(9), 596–602. [DOI: 10.1136/jmg.2008.058628](https://doi.org/10.1136/jmg.2008.058628), [PubMed: 18524835](https://pubmed.ncbi.nlm.nih.gov/18524835/).
5. SanGiovanni J. P., Arking D. E., Iyengar S. K., Elashoff M., Clemons T. E., Reed G. F., Henning A. K., Sivakumaran T. A., Xu X., DeWan A., Agrón E., Rochtchina E., Sue C. M., Wang J. J., Mitchell P., Hoh J., Francis P. J., Klein M. L., Chew E. Y., Chakravarti A. Mitochondrial DNA variants of respiratory complex I that uniquely characterize haplogroup T2 are associated with increased risk of age-related macular degeneration *PLoS One* **2009**, 4(5), e5508. [DOI: 10.1371/journal.pone.0005508](https://doi.org/10.1371/journal.pone.0005508), [PubMed: 19434233](https://pubmed.ncbi.nlm.nih.gov/19434233/).
6. Kussmaul L., Hirst J. The mechanism of superoxide production by NADH:Ubiquinone oxidoreductase (complex I) from bovine heart mitochondria *Proc Natl Acad Sci U S A* **2006**, 103(20), 7607–7612. [DOI: 10.1073/pnas.0510977103](https://doi.org/10.1073/pnas.0510977103), [PubMed: 16682634](https://pubmed.ncbi.nlm.nih.gov/16682634/).